

# Development and Testing of Ventricular Fibrillation Detection Software Module

Artem Pogorelov, Yulia Zavyalova, Alexander Borodin

Petrozavodsk State University  
Department of Computer Science



This project is supported by grant KA179 of Karelia ENPI - joint program of the European Union, Russian Federation and the Republic of Finland



14<sup>th</sup> FRUCT conference  
November 11–15, Helsinki, Finland




# Cardiovascular diseases


- Major health problem
- Greatest risk group
- Constant monitoring




- Smartphone
- Application
- User-friendly
- Customizable ECG viewer
- Analyzer
- Emergency call



# How it works?




# Sudden death





# Requirements for the algorithm



- Real-time
- Reliability
- High speed
- Simplicity



# Phase Space Reconstruction



# Phase Space Reconstruction



# Arrhythmia detection

$$d = \frac{\text{Number of visited boxes}}{\text{Number of all boxes}}$$


If  $d > d_0$ , then classify ECG episode as VF,

$$d_0 = 0.15$$



# Testing of the algorithm

- MIT-BIH
- PhysioNet
- Open-source software
- WFDB-library
- WFDB-applications
- ECG records



# Testing: Normal ECG record

The screenshot shows an IDE window titled "Application Output" with a sub-window "TestingSystem". The output text is as follows:

```
count = 90  
d = 0.05625  
NORMAL  
interval # 188  
count = 99  
d = 0.061875  
NORMAL  
interval # 189  
count = 88  
d = 0.055  
NORMAL  
interval # 190  
count = 67  
d = 0.041875  
NORMAL
```

Overlaid on this is a "TestingSystem" dialog box titled "Select Dir" with a list of files:

- /TestingSystem/records/200.txt
- /TestingSystem/records/207\_V1.txt
- /TestingSystem/records/207\_MLII.txt
- /TestingSystem/records/203.txt
- /TestingSystem/records/114\_V5.txt**

The dialog also has a "Check record!" button. At the bottom of the IDE, a tab bar shows: 1 Issues, 2 Search Results, 3 Application Output, 4 Compile Output, 5 QML/JS Console.



# Testing: VF ECG record

The screenshot shows an IDE window titled "Application Output" with a sub-window "TestingSystem". The output text is as follows:

```
count = 288
d = 0.18
WARNING! Ventricular fibrillation!
interval # 222
count = 256
d = 0.16
WARNING! Ventricular fibrillation!
interval # 223
count = 296
d = 0.185
WARNING! Ventricular fibrillation!
interval # 224
count = 296
d = 0.185
WARNING! Ventricular fibrillation!
```

Overlaid on the IDE is a "TestingSystem" dialog box titled "Select Dir". It contains a list of files with the following paths:

- TestingSystem/records/222\_MLII.txt
- TestingSystem/records/210\_MLII.txt
- TestingSystem/records/210\_V1.txt (highlighted in orange)
- TestingSystem/records/samples106.txt
- TestingSystem/records/222\_V1.txt
- TestingSystem/records/200.txt

The dialog box has a "Check record!" button at the bottom.

At the bottom of the IDE, a tab bar shows the following tabs: Issues, Search Results, Application Output, Compile Output, and QML/JS Console.



# Thanks for your attention!

Try CardiaCare now.

Get this App from Nokia Store



NOKIA

